



BTXpress Cytofusioin®

ELECTROFUSION MEDIUM C TECHNICAL SPECIFICATIONS

FEATURES & BENEFITS

- Optimized for eukaryotic electrofusion applications
- Maximum fusion efficiency, high cell viability
- Stable environment for cell alignment
- Low conductivity means minimal heating
- Physiological pH and Balanced Osmolarity
- Contains no animal products

ADVANCED BUFFER FOR HIGH PERFORMANCE CELL FUSION

BTXpress Cytofusioin® Medium C is an advanced electrofusion buffer designed for use with the BTX Hybrimmune System for high performance cell fusion applications. The low conductivity buffer is specially-formulated to minimize cell turbulence during cell alignment and heating during electrofusion for robust cell fusion efficiency and high cell viability.

BTXpress Cytofusioin® Medium C is sterile filtered from the highest quality non-animal, medical-grade reagents. It is the buffer of choice for many commercial biotech and pharmaceutical companies in their hybridoma generation process for monoclonal antibody discovery.

Алматы (727)345-47-04
Ангарск (395)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

BTXpress Cytofus[®]ion TECHNICAL SPECIFICATIONS

CYTOFUSION[®] MEDIUM SPECIFICATIONS

Volume	500 ml
Osmolarity	270-290 mOsm/L
Conductivity @ 25°C	0.080 ± 0.005 mS/cm
pH	7.2 ± 0.2
Endotoxin	< 0.25 EU/ml
Sterility	Sterile filtered
Storage	2 to 8 °C
Shelf-life	18 months from production date (shipped within 6 months of production)

RECOMMENDATIONS FOR USE

Maintain a sterile environment. Standard aseptic techniques are recommended to avoid contamination during use.

Thorough, repeated wash-in. BTXpress Cytofus[®]ion Medium C, is a low conductivity medium designed for efficient electrofusion. Trace amounts of high conductivity solutions such as PBS or tissue culture growth medium can disrupt the fusion process. Therefore, it is critical to wash the cells thoroughly with BTXpress Cytofus[®]ion Medium C prior to the fusion process. For up to 5 x 10⁷ cells, at least two washes in BTXpress Cytofus[®]ion Medium C are recommended. For more than 5 x 10⁷ cells, at least three washes are recommended.

Thoroughly clean electrofusion chamber. To avoid other sources of ionic contamination, clean the electrofusion chamber after each use and rinse thoroughly with sterile, deionized water.

Room Temperature Electrofusion. For maximum efficiency cell fusion, use BTXpress Cytofus[®]ion Medium C at room temperature. Cell washes prior to the final wash may be carried out at 4°C.

Minimize time in buffer. While Cytofus[®]ion Medium C is non-toxic, it does not contain nutrients to support cell viability over long periods of time. For best results, minimize the time that cells are

suspended in BTXpress Cytofus[®]ion Medium C. It is not recommended that cells remain in BTXpress Cytofus[®]ion Medium C longer than one hour subsequent to the final wash.

5:1 Direct Dilution. Post-electrofusion, Cells in BTXpress Cytofus[®]ion Medium C can be diluted in cell culture medium without washing the cells. A minimum dilution of five parts complete culture medium to one part BTXpress Cytofus[®]ion Medium C is recommended. Alternatively, cells may be washed in growth medium to completely remove BTXpress Cytofus[®]ion Medium C prior to culturing.

WARNINGS & DISCLAIMERS

Do not use if tamper-proof seal is missing or bottle is damaged. Damage to the bottle or deliberate tampering may result in contamination of this product. Check product for clarity before use.

BTXpress Cytofus[®]ion Medium is intended for research and investigational purposes only. It is not intended for human use. This product is not considered to be hazardous based on evaluations made under OSHA Hazard Communication Standard 29 CFR 1910.1200.

ORDERING INFORMATION

Catalog #	Description
47-0001	BTXpress Cytofus [®] ion Electroporation Medium C

STORAGE INFORMATION

Store at 2-8°C after opening. Short term storage (i.e. for shipping) at -20°C to +50°C for 7 days is acceptable. Contents may separate upon freezing. If frozen, mix well before use.

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Тамбов (4752)50-40-97
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +375(257)127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл. почта: bxt@nt-rt.ru || сайт: <https://btx.nt-rt.ru/>